



# Manifolds

- Shark Fin Handle
- Low Torque
- Extended Port Length
- Wide Variety

Customized handles and tubing sets available.

May not be available for all OEM applications.

Single use. Sold packaged sterile. For non-sterile and other packaging options, please contact a Merit OEM representative.

SEE RELATED PRODUCTS



Stopcocks  
253



Tubing  
223









Syringes  
195

## Manifolds with Right Port Orientation and Rotating Luer

							
<b>200 PSI (14 BAR)</b>	2 Port Half Body	3 Port Half Body	4 Port Half Body	4 Port Block Body	3 Port Block Body	2 Port Block Body	1 Port Block Body
ON Handles	<b>202HN-R</b>	<b>203HN-R</b>	<b>204HN-R</b>	<b>204BN-R</b>	<b>203BN-R</b>	<b>202BN-R</b>	<b>201BF-R</b>
OFF Handles	<b>202HF-R</b>	<b>203HF-R</b>	<b>204HF-R</b>	<b>204BF-R</b>	<b>203BF-R</b>	<b>202BF-R</b>	
360/ON Handles					<b>203BN-RT</b>	<b>202BN-RT</b>	
360/OFF Handles			<b>204HF-RT</b>	<b>204BF-RT</b>	<b>203BF-RT</b>		
<b>500 PSI (34 BAR)</b>	2 Port Half Body	3 Port Half Body	4 Port Half Body	4 Port Block Body	3 Port Block Body	2 Port Block Body	1 Port Block Body
ON Handles	<b>502HN-R</b>	<b>503HN-R</b>	<b>504HN-R</b>	<b>504BN-R</b>	<b>503BN-R</b>	<b>502BN-R</b>	<b>501BN-R</b>
OFF Handles	<b>502HF-R</b>	<b>503HF-R</b>	<b>504HF-R</b>	<b>504BF-R</b>	<b>503BF-R</b>	<b>502BF-R</b>	<b>501BF-R</b>
360/ON Handles			<b>202BN-F</b>	<b>504BN-RT</b>			
360/OFF Handles	<b>502HF-RT</b>	<b>503HF-RT</b>		<b>504BF-RT</b>	<b>503BF-RT</b>	<b>502BF-RT</b>	

## Manifolds with Left Port Orientation and Rotating Luer


						
<b>200 PSI (14 BAR)</b>	2 Port Half Body	3 Port Half Body	4 Port Half Body	4 Port Block Body	3 Port Block Body	2 Port Block Body
ON Handles	<b>202HN-L</b>				<b>203BN-L</b>	<b>202BN-L</b>
OFF Handles	<b>202HF-L</b>	<b>203HF-L</b>	<b>204HF-L</b>		<b>203BF-L</b>	<b>202BF-L</b>
<b>500 PSI (34 BAR)</b>	2 Port Half Body	3 Port Half Body	4 Port Half Body	4 Port Block Body	3 Port Block Body	2 Port Block Body
OFF Handles	<b>502HF-L</b>	<b>503HF-L</b>		<b>504BF-L</b>	<b>503BF-L</b>	<b>502BF-L</b>
360/ON Handles	<b>502HN-L</b>	<b>503HN-LT</b>		<b>504BN-LT</b>		
360/OFF Handles	<b>502HF-LT</b>					

# Manifolds with All Female Ports



200 PSI, 5 Port, Block Body,  
360/OFF Red, White & Blue Handles

**500237001**



## 200 PSI (14 BAR)

2 Port Half Body	3 Port Half Body	4 Port Half Body	2 Port Block Body	3 Port Block Body	4 Port Block Body	5 Port Block Body
---------------------	---------------------	---------------------	----------------------	----------------------	----------------------	----------------------

ON Handles

OFF Handles

360/ON Handles

360/OFF Handles

**202HF-F**

**203HF-F**

**202BN-F**

**203BN-F**

**202BF-F**

**205BF-F**

**202HN-FT**

**203HN-FT**

**204HF-FT**

**202BF-FT**

**203BF-FT**

**204BF-FT**

**205BN-FT**

**205BF-FT**

## 500 PSI (34 BAR)

2 Port Half Body	3 Port Half Body	4 Port Half Body	2 Port Block Body	3 Port Block Body	4 Port Block Body	5 Port Block Body
---------------------	---------------------	---------------------	----------------------	----------------------	----------------------	----------------------

ON Handles

OFF Handles

360/ON Handles


**504HN-F**

**503BF-F**

**503BN-FT**

# Manifolds with Right Port Orientation

6" (15cm) tubing extension and rotating luer.



## 200 PSI (14 BAR)

2 Port Half Body	3 Port Half Body	4 Port Half Body	3 Port Block Body	2 Port Block Body
---------------------	---------------------	---------------------	----------------------	----------------------

ON Handles

OFF Handles

360/OFF Handles

**202HN-RX**

**203HN-RX**

**202HF-RX**

**203HF-RX**

**204HF-RX**

**203BF-RX**

**203BF-RXT**

## 500 PSI (34 BAR)

2 Port Half Body	3 Port Half Body	4 Port Half Body	3 Port Block Body	2 Port Block Body
---------------------	---------------------	---------------------	----------------------	----------------------

ON Handles

OFF Handles

360/ON Handles

**503HN-RX**

**502HF-RX**

**503HF-RX**

**504HF-RX**


**503BF-RX**

**502BF-RX**

**504BF-RXT**

## Manifolds with Left Port Orientation

6" (15cm) tubing extension and rotating luer.



**500 PSI (34 BAR)**

3 Port  
Block Body


4 Port  
Block Body

ON Handles  
OFF Handles

**503BN-LX**

**504BF-LX**

## Devos™ Manifolds with Check Relief Valve, Right Port Orientation and Rotating Luer



**200 PSI (14 BAR)**

3 Port  
Half Body

3 Port  
Half Body

2 Port  
Half Body

2 Port  
Half Body

OFF Handles  
360/OFF Handles

**2V3HF-R**

**2V3HF-RT**

**2V2HF-R**

**2V2HF-RT**

Custom colored handles and tubing options available.

