




PERIPHERAL INTERVENTION

Performa[®]

Diagnostic Peripheral Catheters

Push it to the Limit[™]

With Merit's Diagnostic
Peripheral Catheters



**TORQUE
PUSH
FLOW**

Performa[®]

Diagnostic Peripheral Catheters

Introducing Merit Medical's line of Diagnostic Peripheral Catheters.

Designed to provide superior torque, increased shaft strength for better pushability and large inner lumens for improved flow rates.

Now you can **PUSH IT TO THE LIMIT** in today's endovascular environment with Merit's exceptional products.

INNER-LUMEN SIZES

- 4F .042" (1.07mm)
- 5F .046" (1.17mm)
- 6F .054" (1.37mm)
- Yields higher flow rates at lower pressures
- Superior guide wire trackability

STRAIN RELIEF


- Increased strength at catheter/hub interface
- Color coded for easy identification

PRESSURE

- High pressure limits at given flow rates

WIRE BRAID DESIGN

- Greater torque for maximum control
- Improved longitudinal stability
- High degree of maneuverability
- Provides superior pushability



4F CATHETER LINE

- Accepts .038" (0.97Mm) guide wire for physician convenience
- Large inner diameter .042" (1.07Mm) provides excellent flow rates
- Bumper tip added for softness and patient safety

WINGED POLYCARBONATE HUB

- Enhanced handling and control

NYLON PEBAX SHAFT


- Kink Resistant
- Optimal surface smoothness
- Thermal stability at body temperature
- Thin wall design for improved flow rate

RADIOPACITY

- Good shaft visualization under fluoroscopy
- Tip opacity ensures accurate placement and visualization

TIP

- Soft tip for smooth atraumatic positioning
- Excellent tip memory



Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

HEADHUNTER

7715-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7715-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7715-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-50	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7715-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7715-60	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5575-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5575-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	


1



3



5




6



NEWTON

7722-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-40	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7722-E0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7722-50	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5577-E3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5577-53	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	

1



2



3



4



5



BENTSON

7721-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7721-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5578-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5578-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7721-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7721-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5578-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5578-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7721-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7721-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5578-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5578-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	


1



2



3



Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

MANI

7716-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7716-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5579-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5579-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



MANI

VERTEBRAL

7717-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7717-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7717-C0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7717-30	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
7717-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7717-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5583-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5583-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5583-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5583-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



VERTEBRAL

MODIFIED CEREBRAL

7723-A0	4	.042"/1.07mm	90cm	.035" (0.89mm)	-	
7723-10	4	.042"/1.07mm	90cm	.038" (0.97mm)	-	bumper tip
5588-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5588-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



MODIFIED CEREBRAL

BERENSTEIN

7718-C0	4	.042"/1.07mm	40cm	.035" (0.89mm)	-	
7718-30	4	.042"/1.07mm	40cm	.038" (0.97mm)	-	bumper tip
7718-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7718-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7718-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7718-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5592-C3	5	.046"/1.17mm	40cm	.035" (0.89mm)	-	
5592-33	5	.046"/1.17mm	40cm	.038" (0.97mm)	-	
5592-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5592-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5592-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5592-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



BERENSTEIN

STRAIGHT SELECTIVE

7748-A0	4	.042"/1.07mm	110cm	.035" (0.89mm)	-	
7748-13	4	.042"/1.07mm	110cm	.038" (0.97mm)	-	bumper tip




STRAIGHT SELECTIVE

Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

MW2

5590-13	5	.046"/1.17mm	40cm	.038" (0.97mm)	-	
5590-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5590-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	



MODIFIED MW2

7745-A0	4	.042"/1.07mm	40cm	.035" (0.89mm)	2	
7745-10	4	.042"/1.07mm	40cm	.038" (0.97mm)	2	bumper tip
7745-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7745-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	2	



OSBORN

7746-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7746-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
7746-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7746-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7746-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7746-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip



HOOK 0.8

7711-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7711-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5585-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5585-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



HOOK 1.0

7711-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7711-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5585-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5585-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



MODIFIED HOOK

5585-C3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5585-D3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5585-D3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	

1
2
3



Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

SIMMONS

7720-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7720-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5576-F3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5576-A3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5576-13	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7720-B0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7720-20	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5576-G3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5576-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5576-23	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7720-C0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7720-30	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5576-C3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5576-33	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	
7720-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7720-40	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5576-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5576-43	5	.046"/1.17mm	100cm	.038" (0.97mm)	-	

MODIFIED SIMMONS 1.5


7720-F0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
5576-E3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	

MODIFIED SIMMONS

7720-E0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	bumper tip
7720-50	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
7720-70	4	.042"/1.07mm	110cm	.038" (0.97mm)	-	bumper tip
7720-80	4	.042"/1.07mm	120cm	.038" (0.97mm)	-	bumper tip

COBRA 1

7710-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7710-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7710-D0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7710-40	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5584-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5584-D3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5584-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5584-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	
7710-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7710-E0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7710-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7710-50	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5584-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5584-E3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5584-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5584-53	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	
7710-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7710-F0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7710-30	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7710-60	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5584-C3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5584-F3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5584-33	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5584-63	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

SHEPHERD HOOK 0.8

7712-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7712-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5586-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5586-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



SHEPHERD HOOK 1.0

7712-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7712-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
5586-B3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5586-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	




SHEPHERD HOOK 1.0 MODIFIED

7712-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7712-D0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7712-30	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7712-40	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5586-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



RENAL DOUBLE CURVE

7713-A0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7713-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7713-10	4	.042"/1.07mm	65cm	.038" (0.97mm)	2	bumper tip
7713-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
5587-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5587-13	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5587-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	2	



AMIR MOTARJEME®

7747-C0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7747-D0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7747-E0	4	.042"/1.07mm	65cm	.035" (0.89mm)	2	
7747-F0	4	.042"/1.07mm	100cm	.035" (0.89mm)	2	
5591-A3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5591-B3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	
5591-C3	5	.046"/1.17mm	65cm	.035" (0.89mm)	2	
5591-D3	5	.046"/1.17mm	100cm	.035" (0.89mm)	2	



AMIR MOTARJEME CANE

7747-B0	4	.042"/1.07mm	65cm	.035" (0.89mm)	-	
7747-20	4	.042"/1.07mm	65cm	.038" (0.97mm)	-	bumper tip
7747-A0	4	.042"/1.07mm	100cm	.035" (0.89mm)	-	
7747-10	4	.042"/1.07mm	100cm	.038" (0.97mm)	-	bumper tip
5591-E3	5	.046"/1.17mm	65cm	.035" (0.89mm)	-	
5591-F3	5	.046"/1.17mm	100cm	.035" (0.89mm)	-	



Catalog number	French Size	Internal Diameter	Length	Tip Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------	------------------------	-----------------	-------------

HOCKEY STICK

7752-10	4	.042"/1.07mm	100cm	4.5cm	.038" (0.97mm)	-	bumper tip
7543-13	5	.046"/1.17mm	100cm	4.5cm	.038" (0.97mm)	-	
7541-13	6	.054"/1.37mm	100cm	4.5cm	.038" (0.97mm)	-	



HOCKEY STICK

Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

REUTER

7714-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	2	
7714-10	4	.042"/1.07mm	80cm	.038" (0.97mm)	2	bumper tip
7714-B0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7714-20	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
5589-A3	5	.046"/1.17mm	80cm	.035" (0.89mm)	-	
5589-B3	5	.046"/1.17mm	80cm	.035" (0.89mm)	2	
5589-13	5	.046"/1.17mm	80cm	.038" (0.97mm)	-	
5589-23	5	.046"/1.17mm	80cm	.038" (0.97mm)	2	



REUTER

MIKAELSSON

7719-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	-	
7719-10	4	.042"/1.07mm	80cm	.038" (0.97mm)	-	bumper tip
5594-13	5	.046"/1.17mm	80cm	.038" (0.97mm)	-	



MIKAELSSON

KA CATHETERS

5582-23	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
5582-33	5	.052"/1.32mm	65cm	.038" (0.97mm)	-	high flow non-braided



KA CATHETERS

KA2 CATHETERS

5582-13	5	.046"/1.17mm	40cm	.038" (0.97mm)	-	
5582-43	5	.046"/1.17mm	65cm	.038" (0.97mm)	-	
7514-13	6	.054"/1.37mm	40cm	.038" (0.97mm)	-	
7514-23	6	.054"/1.37mm	65cm	.038" (0.97mm)	-	



KA2 CATHETERS

DVS A1 (10cm taper tip)

5595-13	5	.046"/1.12mm	65cm	.035" (0.89mm)	-	
5595-33	5	.046"/1.12mm	100cm	.035" (0.89mm)	-	
5595-23	5	.046"/1.12mm	65cm	.038" (0.97mm)	-	
5595-43	5	.046"/1.12mm	100cm	.038" (0.97mm)	-	



DVS A1

DVS A2 (15cm taper tip)

5595-A3	5	.046"/1.12mm	65cm	.035" (0.89mm)	-	
5595-C3	5	.046"/1.12mm	100cm	.035" (0.89mm)	-	
5595-B3	5	.046"/1.12mm	65cm	.038" (0.97mm)	-	
5595-D3	5	.046"/1.12mm	100cm	.038" (0.97mm)	-	



DVS A2

Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Description
----------------	-------------	-------------------	--------	------------------------	-----------------	-------------

UHF SHEPHERD FLUSH

7602-A1	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	non-braided
7602-B1	5	.052"/1.32mm	90cm	.035" (0.89mm)	6	non-braided
7602-C1	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	non-braided
7602-D1	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	non-braided



ULTRA BOLUS FLUSH

7674-11	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	non-braided
7674-21	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	non-braided
7674-31	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	non-braided
7674-41	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	non-braided




PIGTAIL FLUSH

7749-A0	4	.042"/1.07mm	80cm	.035" (0.89mm)	10	braided
7749-40	4	.042"/1.07mm	80cm	.038" (0.97mm)	10	braided



ULTRA HIGH FLOW PIGTAIL

7602-41	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	non-braided
7602-51	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	non-braided
7602-61	5	.052"/1.32mm	110cm	.035" (0.89mm)	10	non-braided
7602-71	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	non-braided
7602-81	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	non-braided
7602-91	5	.052"/1.32mm	110cm	.038" (0.97mm)	10	non-braided



STRAIGHT FLUSH

7673-11	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	non-braided
7673-21	5	.052"/1.32mm	90cm	.035" (0.89mm)	10	non-braided
7673-31	5	.052"/1.32mm	65cm	.035" (0.89mm)	-	non-braided
7673-41	5	.052"/1.32mm	65cm	.038" (0.97mm)	10	non-braided
7673-51	5	.052"/1.32mm	90cm	.038" (0.97mm)	10	non-braided
7673-61	5	.052"/1.32mm	65cm	.038" (0.97mm)	-	non-braided



MODIFIED HOOK FLUSH

7603-A1	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	non-braided
7603-B1	5	.052"/1.32mm	90cm	.035" (0.89mm)	6	non-braided
7603-C1	5	.052"/1.32mm	65cm	.038" (0.97mm)	6	non-braided
7603-D1	5	.052"/1.32mm	90cm	.038" (0.97mm)	6	non-braided



Catalog number	French Size	Internal Diameter	Length	Recommended Guide wire	# of Side Ports	Bands/Spacing	Tip Style
----------------	-------------	-------------------	--------	------------------------	-----------------	---------------	-----------

PERFORMA VESSEL SIZING All side ports between tip and first band


7602-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 bands 2cm spacing	Pigtail
7602-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 bands 2cm spacing	Pigtail
7674-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 bands 2cm spacing	Ultra Bolus
7674-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 bands 2cm spacing	Ultra Bolus
7673-2M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 bands 2cm spacing	Straight
7673-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 bands 2cm spacing	Straight
7603-2M100	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	2 bands 2cm spacing	Mod. Hook
7603-2M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	2 bands 2cm spacing	Mod. Hook
7742-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Pigtail
7742-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Pigtail
7744-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Straight
7744-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Straight
7743-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 (2cm apart)	Ultra Bolus
7743-2M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	2 (2cm apart)	Ultra Bolus
7755-2M65	4	.041"/1.04mm	65cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook
7755-2M100	4	.041"/1.04mm	100cm	.035" (0.89mm)	6	2 (2cm apart)	Mod. Hook
7602-11M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	11 bands 1cm spacing	Pigtail
7602-11M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	11 bands 1cm spacing	Pigtail
7673-11M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	11 bands 1cm spacing*	Straight
7673-11M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	11 bands 1cm spacing*	Straight
7603-11M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	11 bands 1cm spacing	Mod. Hook
7603-11M100	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	11 bands 1cm spacing	Mod. Hook
7742-20M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	20 bands 1cm spacing	Pigtail
7742-20M100SH	4	.041"/1.04mm	100cm	.035" (0.89mm)	5	20 bands 1cm spacing	Pigtail
7602-20M	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 bands 1cm spacing	Pigtail
7602-20M65	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 bands 1cm spacing	Pigtail
7602-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	5	20 bands 1cm spacing*	Pigtail
7602-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	5	20 bands 1cm spacing*	Pigtail
7673-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 bands 1cm spacing*	Straight
7673-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 bands 1cm spacing*	Straight
7603-20M100SH	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	20 bands 1cm spacing*	Mod. Hook
7603-20M65SH	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	20 bands 1cm spacing*	Mod. Hook



STRAIGHT

PERFORMA VESSEL SIZING All side ports between bands

56535AVS2PIG	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
510035AVS2PIG	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Pigtail
56535AVS2STR	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
510035AVS2STR	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	2 (2cm apart)	Straight
56535AVS2MHK	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	2 (2cm apart)	Mod Hook
510035AVS2MHK	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	2 (2cm apart)	Mod Hook
46535AVS20MHK	4	.041"/1.04mm	65cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook
56535AVS20PIG	5	.052"/1.32mm	65cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
510035AVS20PIG	5	.052"/1.32mm	100cm	.035" (0.89mm)	10	20 (1cm apart)	Pigtail
56535AVS20MHK	5	.052"/1.32mm	65cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook
510035AVS20MHK	5	.052"/1.32mm	100cm	.035" (0.89mm)	6	20 (1cm apart)	Mod Hook




SOFTOUCH PEDIATRIC VESSEL SIZING

4408-2M65	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 bands (1cm apart)	Pigtail
4408-2M80	4	.041"/1.04mm	80cm	.035" (0.89mm)	5	2 bands (1cm apart)	Pigtail
4408-2M65SH	4	.041"/1.04mm	65cm	.035" (0.89mm)	5	2 bands (1cm apart)*	Pigtail
4408-2M80SH	4	.041"/1.04mm	80cm	.035" (0.89mm)	5	2 bands (1cm apart)*	Pigtail

* Bands are 16mm closer to pigtail curve

UNDERSTAND. INNOVATE. DELIVER.

Merit Medical Systems is a worldwide leader in the development, manufacturing, and distribution of quality medical devices used in interventional and diagnostic procedures. With a global work force and cutting-edge manufacturing facilities equipped with the most advanced technology, Merit Medical consistently delivers quality products that are backed by excellent customer service and professional expertise.



US Headquarters, South Jordan, Utah



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.208.4300
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Spain
+34 911238406

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000
(Italiano)
+41 225180035

UK
0800 973 115