

DESIGNED FOR COMFORT & ACCURACY



ReSolve®

LOCKING DRAINAGE CATHETER

ReSolve®

LOCKING DRAINAGE CATHETER

Resolve your toughest drainage procedures with Merit's superior locking drainage catheter. Designed for abscess, nephrostomy and miscellaneous drainage procedures, the ReSolve locking drainage catheter comes packaged with a metal stiffening cannula, a plastic stiffening cannula and a trocar stylet.

Large drainage holes to maximize drainage.

Tapered tip and hydrophilic coating (distal 20cm) for a smooth introduction and less patient discomfort.

Color coded hubs that match standard French sizing charts for easy identification.

Flexible polyurethane material to prevent kinking and ensure excellent catheter rebound.

Pigtail straightener allows for easier insertion of stiffening cannulas

Repositioning Tool (included) for unlocking and repositioning the catheter at any time.

Highly flexible radiopaque marker band, integral to the catheter shaft, forged from platinum.

Catalog Number	French Size	Working Length	Max Guide Wire	Hub Color	Drainage Holes
----------------	-------------	----------------	----------------	-----------	----------------

Non-Marker Band

RLC-6-038	6.5F	25cm	0.038" (0.97mm)	Green	7
RLC-8-038	8.5F	25cm	0.038" (0.97mm)	Blue	8
RLC-10-038	10F	25cm	0.038" (0.97mm)	Purple	8
RLC-12-038	12F	25cm	0.038" (0.97mm)	Yellow	8
RLC-14-038	14F	25cm	0.038" (0.97mm)	Orange	9

Marker Band

RLC-6-038MB	6.5F	25cm	0.038" (0.97mm)	Green	7
RLC-8-038MB	8.5F	25cm	0.038" (0.97mm)	Blue	8
RLC-10-038MB	10F	25cm	0.038" (0.97mm)	Purple	8
RLC-12-038MB	12F	25cm	0.038" (0.97mm)	Yellow	8
RLC-14-038MB	14F	25cm	0.038" (0.97mm)	Orange	9

Packaged 5 per box

RC-RTOOL Repositioning Tool
Packaged 15 per box



ReSolve®

DRAINAGE CATHETER TRAYS

Conveniently containing the ReSolve drainage catheter, MAK-NV, drainage bag, and dilators, the ReSolve trays offer all you need from stick to fixation.

Catalog Number	French Size	20cm Dilator French Size	Marker Band	Securement Device
RLC-NV-6	6.5F	-	no	-
RLC-NV-8	8.5F	8F	no	-
RLC-NV-10	10F	8, 10F	no	-
RLC-NV-12	12F	8, 10, 12F	no	-
RLC-NV-14	14F	8, 10, 12, 14F	no	-
RLCMB-NV-6	6.5F	-	yes	-
RLCMB-NV-8	8.5F	8F	yes	-
RLCMB-NV-10	10F	8, 10F	yes	-
RLCMB-NV-12	12F	8, 10, 12F	yes	-
RLCMB-NV-14	14F	8, 10, 12, 14F	yes	-
RLC-NV-6-SFX	6.5F	-	no	StayFIX® (small)
RLC-NV-8-SFX	8.5F	8F	no	StayFIX (small)
RLC-NV-10-SFX	10F	8, 10F	no	StayFIX (small)
RLC-NV-12-SFX	12F	8, 10, 12F	no	StayFIX (small)
RLC-NV-14-SFX	14F	8, 10, 12, 14F	no	StayFIX (large)
RLCMB-NV-6-SFX	6.5F	-	yes	StayFIX (small)
RLCMB-NV-8-SFX	8.5F	8F	yes	StayFIX (small)
RLCMB-NV-10-SFX	10F	8, 10F	yes	StayFIX (small)
RLCMB-NV-12-SFX	12F	8, 10, 12F	yes	StayFIX (small)
RLCMB-NV-14-SFX	14F	8, 10, 12, 14F	yes	StayFIX (large)
RLC-NV-6-REV	6.5F	-	no	Revolution™
RLC-NV-8-REV	8.5F	8F	no	Revolution
RLC-NV-10-REV	10F	8, 10F	no	Revolution
RLC-NV-12-REV	12F	8, 10, 12F	no	Revolution
RLC-NV-14-REV	14F	8, 10, 12, 14F	no	Revolution
RLCMB-NV-6-REV	6.5F	-	yes	Revolution
RLCMB-NV-8-REV	8.5F	8F	yes	Revolution
RLCMB-NV-10-REV	10F	8, 10F	yes	Revolution
RLCMB-NV-12-REV	12F	8, 10, 12F	yes	Revolution
RLCMB-NV-14-REV	14F	8, 10, 12, 14F	yes	Revolution

Packaged as individual units

All kits contain ReSolve Locking Drainage Catheter, MAK-NV-001, and FZ624. Items in **BOLD** are ready for immediate delivery. Allow 6-8 weeks for all others.





LOCKING DRAINAGE CATHETER WITH SECUREMENT DEVICE

Having your drainage catheter and securement device conveniently packaged together saves both time and inventory requirements.

Catalog Number	French Size	Working Length	Marker Band	Securement Device
Non-Marker Band				
RLC-6-REV	6.5F	25cm	no	Revolution™
RLC-6-SFX	6.5F	25cm	no	StayFIX® (small)
RLC-8-REV	8.5F	25cm	no	Revolution
RLC-8-SFX	8.5F	25cm	no	StayFIX (small)
RLC-10-REV	10F	25cm	no	Revolution
RLC-10-SFX	10F	25cm	no	StayFIX (small)
RLC-12-REV	12F	25cm	no	Revolution
RLC-12-SFX	12F	25cm	no	StayFIX (small)
RLC-14-REV	14F	25cm	no	Revolution
RLC-14-SFX	14F	25cm	no	StayFIX (large)
Marker Band				
RLCMB-6-REV	6.5F	25cm	yes	Revolution
RLCMB-6-SFX	6.5F	25cm	yes	StayFIX (small)
RLCMB-8-REV	8.5F	25cm	yes	Revolution
RLCMB-8-SFX	8.5F	25cm	yes	StayFIX (small)
RLCMB-10-REV	10F	25cm	yes	Revolution
RLCMB-10-SFX	10F	25cm	yes	StayFIX (small)
RLCMB-12-REV	12F	25cm	yes	Revolution
RLCMB-12-SFX	12F	25cm	yes	StayFIX (small)
RLCMB-14-REV	14F	25cm	yes	Revolution
RLCMB-14-SFX	14F	25cm	yes	StayFIX (large)

RC-RTOOL Repositioning Tool
Packaged 15 per box



ReSolve®

NON-LOCKING DRAINAGE CATHETER

Resolve your toughest temporary drainage procedures with Merit Medical's non-locking drainage catheter. Designed for temporary drainage of abscess, nephrostomy and miscellaneous drainage procedures, the ReSolve comes in a J or straight tip and is packaged with a metal stiffening cannula and a trocar stylet.

Large drainage holes to maximize drainage.

Tapered tip and hydrophilic coating (distal 10cm) for a smooth introduction and less patient discomfort.

Color-coded hubs that match standard French sizing charts for easy identification.

Flexible polyurethane material to prevent kinking and ensure excellent catheter rebound.



Catalog Number	French Size	Length	Max Guide Wire	Tip Config.	Drainage Holes
RNL-6-038S	6.5F	20cm	0.038" (0.97mm)	Straight	4
RNL-8-038S	8.5F	25cm	0.038" (0.97mm)	Straight	4
RNL-10-038S	10F	25cm	0.038" (0.97mm)	Straight	4
RNL-12-038S	12F	25cm	0.038" (0.97mm)	Straight	5
RNL-6-038J	6.5F	20cm	0.038" (0.97mm)	J Curve	4
RNL-8-038J	8.5F	25cm	0.038" (0.97mm)	J Curve	4
RNL-10-038J	10F	25cm	0.038" (0.97mm)	J Curve	4
RNL-12-038J	12F	25cm	0.038" (0.97mm)	J Curve	5

Packaged 5 per box



DRAINAGE BAGS

These bags come with a VELCRO® brand strap and are designed with patient mobility and comfort in mind.

- Adjustable length tubing to fit a variety of patients
- Cloth backing for patient comfort
- Twist valve for easy disposal of fluids
- Easy-to-read markings to identify fluid volume
- Internal anti-reflux valve prevents backflow

MDD600 600 mL bag

MDD600PML 600 mL bag

FZ624 600 mL cloth-backed bag

Packaged 20 per box

K08-01907 Jackson-Pratt with 100 mL bulb

K10-03450 MDD600 drainage bag with Jackson-Pratt® 100 mL bulb designed for vacuum drainage

Packaged 10 per box

ABSCESS DRAINAGE KIT

K10-00230B OSHA-compliant closed waste system with 1400 mL MDD bag, tubing and CRV Valve
Designed for quick and safe drainage elimination.

Packaged 15 per box

DRAINAGE CONNECTOR TUBING

K10-04129 30 cm connecting tubing

K10-04130 42 cm connecting tubing

K10-04131 30 cm connecting tubing with one-way stopcock

K10-04132 42 cm connecting tubing with one-way stopcock

Packaged 25 per box

DRAINAGE TUBING

FAS5248 48" large bore tubing with pinch clamp and non-vented spike. Packaged 25 units per box.

K08-02027 Tubing kit with barbed male luer. Packaged 25 units per box.

K10-04750 48" Large bore tubing with roller clamp and 16G x 1.5" (38.1 mm) needle.
Packaged 15 per box

VALVE ADAPTER

VA-40 Designed to be back loaded over a guide wire and attached to the MAK-NV to facilitate aspiration or infusion. Accepts up to a 4F device

Packaged 10 per box

LARGE BORE STOPCOCKS

Clear body with large lumen designed to increase flow rate.

L3SNC 3-way large bore stopcock

L4SNC 4-way large bore stopcock

Packaged 50 per box

VacLok® SYRINGES

- Patented clear body design with multiple locking positions.
- Maintains negative pressure while reducing physician strain.

VAC120 20 mL VacLok Syringe

VAC130 30 mL VacLok Syringe


Packaged 20 per box

VAC110 10 mL VacLok Syringe

VAC160 60 mL VacLok Syringe

Packaged 25 per box

VELCRO® brand strap is a registered trademark of Velcro Industries B.V.
Jackson-Pratt® is a registered trademark of Cardinal Health.



ReSolve®

LOCKING DRAINAGE CATHETER

To prevent irritation of the surrounding tissues, the tip of the ReSolve Locking Drainage Catheter is positioned within the pigtail curve. Depth markers located on the catheter shaft aid accurate placement and act as a reference for correct positioning.



Features and Benefits

Accurate placement

Depth markers aid in accurate placement and act as a reference for correct positioning.

Atraumatic

Catheter tip positioned within pigtail curve to reduce tissue irritation.

Ease of Insertion

Tapered tip and hydrophilic coating facilitate a smooth introduction to reduce patient discomfort.

Catheter Security

The unique design makes it easy to engage the locking mechanism with one hand and is also tamper resistant. The low profile of the locking mechanism and hub are designed to provide patient comfort.

Drainage Efficiency

The large drainage holes maximize drainage while maintaining catheter patency.

High Visibility

Radiopaque platinum marker band incorporated into the catheter wall for visualizing the most proximal drainage hole.

Easy Identification

Color-coded catheter hubs and packaging make it easy to determine the French size for initial placements and catheter exchanges.



Ordering Information

Catalog Number	French Size	Working Length	Marker Band	Hub Color	Drainage Holes
RLC15-6-038MB	6.5F	15 cm	Yes	Green	8
RLC15-8-038MB	8.5F	15 cm	Yes	Blue	8
RLC15-10-038MB	10F	15 cm	Yes	Purple	8
RLC15-12-038MB	12F	15 cm	Yes	Yellow	9
RLC15-14-038MB	14F	15 cm	Yes	Orange	9
RLC-6-038MB	6.5F	25 cm	Yes	Green	8
RLC-8-038MB	8.5F	25 cm	Yes	Blue	8
RLC-10-038MB	10F	25 cm	Yes	Purple	8
RLC-12-038MB	12F	25 cm	Yes	Yellow	9
RLC-14-038MB	14F	25 cm	Yes	Orange	9
RLC40-6-038MB	6.5F	40 cm	Yes	Green	8
RLC40-8-038MB	8.5F	40 cm	Yes	Blue	8
RLC40-10-038MB	10F	40 cm	Yes	Purple	8
RLC40-12-038MB	12F	40 cm	Yes	Yellow	9
RLC40-14-038MB	14F	40 cm	Yes	Orange	9
RLC-6-038	6.5F	25 cm	-	Green	8
RLC-8-038	8.5F	25 cm	-	Blue	8
RLC-10-038	10F	25 cm	-	Purple	8
RLC-12-038	12F	25 cm	-	Yellow	9
RLC-14-038	14F	25 cm	-	Orange	9

Packaged 5 per box

Contact Merit Medical for more information on our line of drainage products.

Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

www.merit.com

Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) • Amerikalaan 42, 6199 AE Maastricht-Airport • The Netherlands • Tel: +31 43 358 82 22

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24 • France 0800 91 60 30

Finland 0800 770 586 • Germany 0800 182 0871 • Ireland (republic) 1800 553 163 • Italy 800 897 005 • Luxembourg 8002 25 22 • Netherlands 0800 022 81 84

Norway 800 11629 • Sweden 020 792 445 • UK 0800 973 115

402815001/A ID 121212